

## Joerns Custom Sling Order Form Full Back

STEP 1 Co	ontact Detai	B				
Company:						<u> </u>
Delivery Address:						
Tel:					. •	
Email:						- 1
			doubt, please co		's standard specifica	ation on page 4.
	□ xs	☐ Small	☐ Medium	☐ Large	☐ XL	Custom
Dimension A	810mm	980mm	1180mm	1360mm	1690mm	mm
Dimension B	860mm	890mm	1060mm	1160mm	1340mm	mm
Dimension C	1310mm	1370mm	1580mm	1760mm	2090mm	mm
Note: Some slings inco	Polyester Solid Polyester Net	one fabric type - pleas	Super-Soft  3D Spacer  se consider this when		Instructions:  al changes. See page 4 fo	or more information.
STEP 4 Pa	ndding Optio	Ons: Add padding	to your sling for extra	comfort.		
	5mm HD Foam		m Soft Foam nm Soft Foam		3mm Spacer	☐ Hollow Fibre
Head Support Body Legs	Special Inst	tructions:				
STEP 5 In	ner Lining C	ptions: A lining	g can be added to the	inside of your sling	g for enhanced comfort.	
	Micro Fleece	Fixed (Stitched)	Detachable (Velcro®)	☐ Hea	ad Body	☐ Legs
Special Instructions:						

STEP 6	Head Sup	port: Use page	3 to help illustrate th	e shape and size of h	nead support if required.		
Adjust Head Support		Special Instructions:					
STEP 7	Straps, Lo	ops & Han	dles: Adjust stra	p lengths, number of	loops, loop colours and positioning handles.		
	Standard Strap Length	Custom Strap Length	Standard No. of Loops	No. of Extra Loops?	Special Instructions: (e.g. quantity and placement of positioning handles)		
Shoulder	440mm	mm	4				
Leg	440mm	mm	2				
☐ Add Posit	ioning Handle(s)		Omit Modesty Lo	ор			
☐ Omit Posi	itioning Handle(s)						
					election/sequence in the special instructions box above.		
See page 3 for dir	mensional illustration	ons for standard str	ap lengths and loop	positions.			
STEP 8	Body Belt	Add or alter a boo	ly belt to enhance s	upport and security d	uring transfer. Use page 3 to help illustrate your needs.		
☐ Buckle	Belt Lengt	n Belt Wid	Special li	nstructions: (e.g. pade	ding, elastic, etc)		
☐ Velcro®	mn	1	mm				
STEP 9	Other Opt	ions: Use page	3 to help illustrate y	our bespoke requiren	nents.		
☐ Add Positi	ioning Pocket	☐ Omit Centr	e Markers	pecial Instructions:			
STEP 10	Additiona	l Requirem	ents				
STEP 11	Confirmat	ion					
☐ I am happ	y with and hereby	give permission to p	proceed with the cus	stom sling order defin	ed on this form.		
☐ Please cal	Il me to discuss/co	nfirm this order and	or price prior to ma	nufacture. Ph	one No:		
Customer PO	Number:				red By:		
	Quantity:			· · · · · · · · · · · · · · · · · · ·	Signed:		
Sales Order	Number:				Date:		

Please Note: Joerns Healthcare reserve the right to refuse a special order if deemed unsafe. In such event, we will contact you to discuss the matter further.

## Please return your completed form to:


orders@joerns.co.uk

Please attach any supporting illustrations when submitting your form.

**(T)** +44 (0)344 811 1158 **(W)** www.joerns.co.uk

## Illustrative Guidance

Please use this page to help illustrate or clarify any bespoke requirements



## Standard Specification Breakdown

Please consider the following information when requesting changes to your standard sling.

	Sling Model					
	Oxford Full Back	Oxford Full Back Net	Oxford Full Back Spacer			
Material						
Core Material	Polyester Solid (Grey)	Polyester Net (Grey)	3D Spacer (Black)			
Spine Overlay	N/A	N/A	N/A			
Leg Overlay	N/A	N/A	N/A			
Padding						
Head Support	None	None	None			
Body	None	None	None			
Legs	5mm HD Foam	5mm HD Foam	5mm Soft Foam			
Head Support						
Included	✓	✓	✓			
Integrated Stiffeners	✓	✓	✓			
Adjustment Straps	N/A	N/A	N/A			
Loops, Clips & Handles						
Sling Attachment Type	Loops	Loops	Loops			
Positioning Handles	2	2	3			
Modesty Loop	✓	✓	✓			
Body Belt						
Included	N/A	N/A	N/A			
Attachment Type	N/A	N/A	N/A			
Other Options						
Centre Markers	Centre Markers ✓		✓			
Positioning Pocket	N/A	N/A	N/A			

You do not need to submit this page, it is for your reference only.