


# Custom Sling Order Form

## High Back

### STEP 1 ➤ Contact Details:

Company:	
Delivery Address:	
Tel:	
Email:	



Prior to requesting changes or if you're in any doubt, please consult this sling's standard specification on page 4.

### STEP 2 ➤ Core Sling Dimensions: Select a standard size or enter custom dimensions.

	<input type="checkbox"/> XS	<input type="checkbox"/> Small	<input type="checkbox"/> Medium	<input type="checkbox"/> Large	<input type="checkbox"/> XL	Custom
Dimension A	720mm	860mm	990mm	1040mm		mm
Dimension B	610mm	750mm	840mm	840mm		mm
Dimension C	890mm	1100mm	1250mm	1320mm		mm

**Note:** Use page 3 to illustrate bespoke sizing. Dimension A is the mid-point of the sling body. All dimensions are subject to manufacturing tolerances.

### STEP 3 ➤ Material Options: Select the core material for your sling and include any bespoke requirements in the box provided.




**Note:** Some slings incorporate more than one fabric type - please consider this when requesting material changes. See page 4 for more information.

### STEP 4 ➤ Padding Options: Add padding to your sling for extra comfort.



### STEP 5 ➤ Inner Lining Options: A lining can be added to the inside of your sling for enhanced comfort.



## STEP 6 ➤ Head Support: Use page 3 to help illustrate the shape and size of head support if required.

Adjust Head Support

Special Instructions:

## STEP 7 ➤ Straps, Loops, Clips & Handles: Adjust strap lengths, number of loops, loop colours and positioning handles.

	Standard Strap Length	Custom Strap Length	Standard No. of Loops	No. of Extra Loops?
Shoulder	280mm	mm	2	
Leg	800mm	mm	2	

Special Instructions: (e.g. quantity and placement of positioning handles)

Add Positioning Handle(s)

Clip Attachments (Securi3®)

When choosing clips, please specify the length of strap required from the sling body to the clip (required for both shoulder and leg attachment points).

Omit Positioning Handle(s)

Add Modesty Loop

**Loop Colour Options:** Black, blue, green, grey, orange, red, white and yellow. Please confirm selection/sequence in the special instructions box above. See page 3 for dimensional illustrations for standard strap lengths and loop positions.

## STEP 8 ➤ Body Belt: Add or alter a body belt to enhance support and security during transfer. Use page 3 to help illustrate your needs.

Buckle

**Belt Length**

**Belt Width**

Special Instructions: (e.g. padding, elastic, etc)

Velcro®

## STEP 9 ➤ Other Options: Use page 3 to help illustrate your bespoke requirements.

Omit Positioning Pocket

Omit Centre Markers

Special Instructions:

## STEP 10 ➤ Additional Requirements

## STEP 11 ➤ Confirmation

I am happy with and hereby give permission to proceed with the custom sling order defined on this form.

Please call me to discuss/confirm this order and/or price prior to manufacture.

Phone No:

Customer PO Number:	
Order Quantity:	
Sales Order Number:	

Ordered By:	
Signed:	
Date:	

**Please Note:** Joerns Healthcare reserve the right to refuse a special order if deemed unsafe. In such event, we will contact you to discuss the matter further.

**Please return your completed form to:**


**orders@joerns.co.uk**

**Please attach any supporting illustrations when submitting your form.**

**(T) +44 (0)344 811 1158 (W) www.joerns.co.uk**

## Illustrative Guidance

Please use this page to help illustrate or clarify any bespoke requirements



# Standard Specification Breakdown

Please consider the following information when requesting changes to your standard sling.

Sling Model	
Oxford High Back	
<b>Material</b>	
Core Material	Polyester Solid (Black)
Spine Overlay	N/A
Leg Overlay	N/A
<b>Padding</b>	
Head Support	None
Body	None
Legs	5mm HD Foam
<b>Head Support</b>	
Included	✓
Integrated Stiffeners	✓
Adjustment Straps	N/A
<b>Loops, Clips &amp; Handles</b>	
Sling Attachment Type	Loops
Positioning Handles	N/A
Modesty Loop	N/A
<b>Body Belt</b>	
Included	N/A
Attachment Type	N/A
<b>Other Options</b>	
Centre Markers	✓
Positioning Pocket	✓

You do not need to submit this page, it is for your reference only.