O	ns ford ®	Custom Sling O Long Seat	rder Form		/
STEP 1 Co	ontact Details:			в	
Company:					
Delivery Address:				<▲	
Tel:					
Email:				7	•

Prior to requesting changes or if you're in any doubt, please consult this sling's standard specification on page 4.

STEP 2 Core SI	ing Dimensions:	Select a standard size or ent	er custom dimensions.			
	Small	Medium	Large	Custom		
Dimension A	590mm	690mm	730mm	mm		
Dimension B	1020mm	1170mm	1280mm	mm		
			dy. All dimensions are subject to n nd include any bespoke requiremen			
	ster Solid	Super-Soft	Special Instructions:			
	ster Net	☐ 3D Spacer				
Note: Some slings incorporate i	Note: Some slings incorporate more than one fabric type - please consider this when requesting material changes. See page 4 for more information.					
STEP 4 Padding		ing to your sling for extra com 5mm Soft Foam 10mm Soft Foam	fort.	Hollow Fibre		
 Head Sup. Section Body Leg Section 	Special Instructions:					
STEP 5 Inner Li Micro Fleec Special Instructions:	Fixed	ning can be added to the insid	le of your sling for enhanced comfo	ort. ody Leg Section		
L						

STEP 6 Head Support: Use page 3 to help illustrate the shape and size of head support if required.

 Adjust Head Support Section
 Special Instructions:

STEP 7 Straps, Loops & Handles: Adjust strap lengths, number of loops, loop colours and positioning handles.

		Standard Strap Length	Custom Strap Length	Standard No. of Loops	No. of Extra Loops?	Special Instructions: (e.g. quantity and placement of positioning handles)
Sho	oulder	210mm	mm	2		
	Hip	640mm	mm	2		
	Leg	310mm	mm	2		
Add Positioning Handle(s)						

Loop Colour Options: Black, blue, green, grey, orange, red, white and yellow. Please confirm selection/sequence in the special instructions box above. See page 3 for dimensional illustrations for standard strap lengths and loop positions.

STE	Body Belt: Add or alter a body belt to enhance support and security during transfer. Use page 3 to help illustrate your needs.					
	Buckle Belt Length Belt Width Special Instructions: (e.g. padding, elastic, etc) Velcro® mm mm					
STE	STEP 9 Other Options: Use page 3 to help illustrate your bespoke requirements.					
	Add Commode Aperture Omit Centre Markers Special Instructions:					
STEP 11 Confirmation						
	I am happy with and hereby give permission to proceed with the custom sling order defined on this form.					
	Please call me to discuss/confirm this order and/or price prior to manufacture. Phone No:					

Customer PO Number:	Ordered By:	
Order Quantity:	Signed:	
Sales Order Number:	Date:	


Please Note: Joerns Healthcare reserve the right to refuse a special order if deemed unsafe. In such event, we will contact you to discuss the matter further.

Please return your completed form to:

Orders@joerns.co.uk Please attach any supporting illustrations when submitting your form. (T) +44 (0)344 811 1158 (W) www.joerns.co.uk

Illustrative Guidance

Please use this page to help illustrate or clarify any bespoke requirements



Standard Specification Breakdown

Please consider the following information when requesting changes to your standard sling.

	Sling Model		
	Oxford Long Seat	Oxford Long Seat Net	
Material			
Core Material	Polyester Solid (Grey)	Polyester Net (Grey)	
Spine Overlay	N/A	N/A	
Leg Overlay	N/A	N/A	
Padding			
Head Support Section	None	None	
Body	None	None	
Leg Section	None	None	
Head Support			
Included	\checkmark	\checkmark	
Integrated Stiffeners	N/A	N/A	
Adjustment Straps	N/A	N/A	
Loops, Clips & Handles			
Sling Attachment Type	Loops	Loops	
Positioning Handles	2	2	
Modesty Loop	N/A	N/A	
Body Belt			
Included	N/A	N/A	
Attachment Type	N/A	N/A	
Other Options			
Centre Markers	\checkmark	\checkmark	
Positioning Pocket	N/A	N/A	

You do not need to submit this page, it is for your reference only.