# Joerns Custom Sling Order Form UltraFine Reflex

STEP 3



#### Prior to requesting changes or if you're in any doubt, please consult this sling's standard specification on page 4.

<b>STEP 2 Core Sling Dimensions:</b> Select a standard size or enter custom dimensions.							
	🗌 XS	Small	Medium	Large	🗌 XL	Custom	
Dimension A	740mm	820mm	920mm	1000mm	1080mm	mm	
Dimension B	730mm	810mm	900mm	970mm	1040mm	mm	
Dimension C	1060mm	1180mm	1330mm	1420mm	1530mm	mm	

Note: Use page 3 to illustrate bespoke sizing. Dimension A is the mid-point of the sling body. All dimensions are subject to manufacturing tolerances.

Material Options: Select the core material for your sling and include any bespoke requirements in the box provided.

Polyester Solid	Super-Soft	Special Instructions:
Polyester Net	3D Spacer	

Note: Some slings incorporate more than one fabric type - please consider this when requesting material changes. See page 4 for more information.



STEP 6 Head Support: Use page 3 to help illustrate the shape and size of head support if required.

Special Instructions:

Adjust Head Support

STEP 7 Straps, Loops, Clips & Handles: Adjust strap lengths, number of loops, loop colours and positioning handles.

	Standard Strap Length	Custom Strap Length	Standard No. of Loops	No. of Extra Loops?	Special Instructions: (e.g. quantity and placement of positioning handles)
Shoulder	600mm	mm	1		
Hip	N/A	mm	N/A		
Leg	590mm	mm	3		
Add Positioning Handle(s)   Omit Positioning Handle(s)   Omit Modesty Loop			Clip Attachments When choosing clip the length of strap r sling body to the clip shoulder and leg att	s, please specify equired from the o (required for both	

Loop Colour Options: Black, blue, green, grey, orange, red, white and yellow. Please confirm selection/sequence in the special instructions box above. See page 3 for dimensional illustrations for standard strap lengths and loop positions.

<b>STEP 8 Body Belt:</b> Add or alter a body belt to	enhance support and security during transfer. Use page 3 to help illustrate your needs.				
Buckle Belt Length Belt Width	Special Instructions: (e.g. padding, elastic, etc)				
Velcro®   mm   mm					
STEP 9 Other Options: Use page 3 to help   Add Positioning Pocket Add Centre Markers   STEP 10 Additional Requirements	Questicities				

### STEP 11 Confirmation

	I am happy with and hereby give permission to proceed with the custom sling order defined on this form.				
Please call me to discuss/confirm this order and/or price prior to manufacture. Phone No:					

Customer PO Number:	Ordered By:	
Order Quantity:	Signed:	
Sales Order Number:	Date:	

Please Note: Joerns Healthcare reserve the right to refuse a special order if deemed unsafe. In such event, we will contact you to discuss the matter further.

#### Please return your completed form to:

Orders@joerns.co.uk Please attach any supporting illustrations when submitting your form. (T) +44 (0)344 811 1158 (W) www.joerns.co.uk



## Standard Specification Breakdown

Please consider the following information when requesting changes to your standard sling.

	Sling Model
	Oxford UltraFine Reflex
Material	
Core Material	3D Spacer (Black)
Spine Overlay	Rip-stop Nylon (Black)
Leg Overlay	Rip-stop Nylon (Black)
Padding	
Head Support	Hollow Fibre
Body	3D Spacer (Black) to Spine
Legs	3D Spacer (Black)
Head Support	
Included	$\checkmark$
Integrated Stiffeners	✓
Adjustment Straps	N/A
Loops, Clips & Handles	
Sling Attachment Type	Loops
Positioning Handles	1
Modesty Loop	$\checkmark$
Body Belt	
Included	$\checkmark$
Attachment Type	Buckle
Other Options	
Centre Markers	N/A
Positioning Pocket	N/A

You do not need to submit this page, it is for your reference only.