

## **Custom Sling Order Form** UltraFine Straight Leg

(Including Lite)

STEP 1	<b>Contact Detai</b>	ls:	4A			
Company:					C	
Delivery Addres	s: 				-	7 \ \ \
						/
Tel:						1
Email:						1
Prior to r	equesting changes	or if you're in any	doubt, please co	nsult this sling's	s standard speci	fication on page 4.
		_				
STEP 2	Core Sling Di	mensions: Se	elect a standard size o	enter custom dime	ensions.	
	□ XS	☐ Small	☐ Medium	☐ Large	☐ XL	Custom
Dimension		880mm	970mm	1070mm	1180mm	mm
Dimension		610mm	650mm	710mm	800mm	mm
Dimension	950mm	1080mm	1170mm	1290mm	1430mm	mm
	☐ Polyester Solid		☐ Super-Soft	Special I	nstructions:	
	Polyester Net		☐ 3D Spacer			
Note: Some slings	s incorporate more than o	one fabric type - plea	se consider this when	requesting materia	l changes. See page	4 for more information.
STEP 4	Padding Option	ons: Add padding	to your sling for extra	comfort.		
	5mm HD Foam		m Soft Foam		3mm Spacer	☐ Hollow Fibre
☐ Body	Special Ins	tructions:				
Legs						
_						
STEP 5	Inner Lining C	Options: A lining	g can be added to the	nside of your sling	for enhanced comfo	rt.
	☐ Micro Fleece	Fixed (Stitched)	□ Detachable (Velcro®)	е	Body	☐ Legs
Special Instruction	ns:					

This sling is available with head support as standard. If you would like a head support included, please use form CSOF00028 for the UltraFine Straight Leg with Head Support. For a copy of this form, please contact Joerns Healthcare.									
	Otropo Lo	ana Oliv		0 110	_				
STEP 7	Straps, Lo	ops, Cii	os	& Handles	Adjust strap le		er of loops, loop colours and positioning handles.		
	Standard Strap Length			Standard No. of Loops	No. of Extra Loops?		Special Instructions: (e.g. quantity and placement of positioning handles)		
Shoulder	300mm	mm		4					
Hip	N/A	mm		N/A					
Leg	550mm	mr	n	4					
Add Pos	Add Positioning Handle(s)			Clip Attachment	(Securi3®)				
☐ Omit Po	Omit Positioning Handle(s)			When choosing cl					
Omit Mo	☐ Omit Modesty Loop			sling body to the clip (required for both shoulder and leg attachment points).					
See page 3 for o	dimensional illustratio	ons for standar	d str	ap lengths and loo	p positions.		n/sequence in the special instructions box above.		
STEP 8	STEP 8 Body Belt: Add or alter a body belt to enhance support and security during transfer. Use page 3 to help illustrate your needs.								
Belt Length Belt Width Special II					Instructions: (e.	.g. padding, e	lastic, etc)		
☐ Velcro®	mn	mm mm							
STEP 9	Other Opt	ions: Use p	page	3 to help illustrate	your bespoke re	equirements.			
☐ Add Pos	Add Positioning Pocket Add Centre Markers				Special Inst	tructions:			
STEP 10 Additional Requirements									
STEP 11	➤ Confirmat	ion							
☐ I am happy with and hereby give permission to proceed with the custom sling order defined on this form.									
Please call me to discuss/confirm this order and/or price prior to manufacture.									
	_								
Customer PC	O Number:					Ordered By	:		

STEP 6 Head Support:

Order Quantity:

Sales Order Number:

Please Note: Joerns Healthcare reserve the right to refuse a special order if deemed unsafe. In such event, we will contact you to discuss the matter further.

Signed:

Date:

## Please return your completed form to:

orders@joerns.co.uk

Please attach any supporting illustrations when submitting your form.

**(T)** +44 (0)344 811 1158 **(W)** www.joerns.co.uk

## Illustrative Guidance Please use this page to help illustrate or clarify any bespoke requirements Sconn PADMIM

## Standard Specification Breakdown

Please consider the following information when requesting changes to your standard sling.

	Sling Model					
	Oxford UltraFine Straight Leg	Oxford UltraFine Lite Straight Leg				
Material						
Core Material	3D Spacer (Black)	Super-soft (Blue)				
Spine Overlay	N/A	N/A				
Leg Overlay	Rip-stop Nylon (Black)	Rip-stop Nylon (Black)				
Padding						
Head Support	N/A	N/A				
Body	None	None				
Legs	3D Spacer (Black)	3D Spacer (Black)				
Head Support						
Included	N/A	N/A				
Integrated Stiffeners	N/A	N/A				
Adjustment Straps	N/A	N/A				
Loops, Clips & Handles						
Sling Attachment Type	Loops	Loops				
Positioning Handles	N/A	N/A				
Modesty Loop	✓	✓				
Body Belt						
Included	N/A	N/A				
Attachment Type	N/A	N/A				
Other Options						
Centre Markers	N/A	N/A				
Positioning Pocket	N/A	N/A				

You do not need to submit this page, it is for your reference only.